

ESSENTIAL GUIDE WEBCAST

LC|GC's **CHROM**academy
powered by crawford scientific

Introduction to Ion Chromatography

Sponsored by

Thermo
SCIENTIFIC

LIVE WEBCAST

Thursday 29th Nov. 2012, 11:00 AM EST; 15:00 GMT; 16:00 BST

Webcast Notes

- Type your questions in the “Submit Question” box, located below the slide window
- You can enlarge the slide window at any time by clicking on the “Enlarge Slides” button, located below the presentation window – the slides will advance automatically throughout the event
- If you are experiencing technical problems with viewing or hearing the event, please click on the “Help” button, located below the slide window

Speakers



Tony Taylor
Technical Director
Crawford Scientific



Dr Ken Cook
Bio-separation Manager
ThermoFisher Scientific

Moderator



Alasdair Matheson
Editor
LCGC Europe



The Essential Guide **tutorials**

monthly comprehensive
study guides with interactive
content

The Essential Guide **webcasts**

the CHROMacademy archive
is growing fast and each
month our experts create new
and engaging content

troubleshooting

our interactive tools
and calculators help
solve problems

ask the expert

CHROMacademy members have
direct access to our professional
training team. We promise to help
solve your chromatographic
problems

industry news

keeping you
up to date

articles

from John Dolan,
John Hinshaw,
Michael Balogh,
Ron Majors and
many more...

application notes

we have
1000s of articles and
applications

eLearning

CHROMacademy is packed full of
HPLC, GC, Sample Prep and Mass Spec
content – most subjects covered in full
multimedia. Study online or offline / at
work or at home

eLabs

interact with our
virtual machines

assessments

test your knowledge

Outline

- Basic Theory of Ion Chromatography
- Columns
 - Ion Exchange Process
 - Altering Selectivity with buffers and column choice
- Detection and Suppression
 - Basic concept
 - State of the Art
- Modern Day Ion Chromatography
 - Concept and schematics of new Systems
 - Fast Ion Chromatography
 - Capillary Ion Chromatography
 - IC MS

Difficulties associated with the measurement of small ionic species using liquid chromatography

- No retention on traditional reverse phase columns for separation
- Mostly no UV absorbance for detection
- Ion exchange eluents for separation give too high a background to use conductivity as a means of detection

History of Ion Chromatography (IC)

Early 1970s:










Ion chromatography (IC) was developed by Hamish Small and co-workers at Dow Chemical Company as a novel method of IEC usable in automated analysis.

The Dow Chemical Company technology was acquired by Durrum Instrument Corp., which later formed a separate business unit for its new IC products, naming it **Dionex** (**D**ow **Ion** **E**xchange).

Ion chromatography (IC) is a process that allows the separation of ions and polar molecules based on their charge properties. Detection of the separated ions is usually achieved with conductivity detection.


Ion Chromatography Analytes

Functional groups detected with suppressed conductivity
Suppressed conductivity detection




Who uses Ion Chromatography (IC)?


- Environmental Testing Labs
- Water providers
- Wastewater treatment facilities
- Pharmaceutical companies
- Chemical/petrochemical companies
- Food producers and processors
- Power utilities
- Electronics manufacturers
- Mining/metals/plating companies
- Government agencies
- Universities
- and more...




How does Ion Chromatography work?



The Ion Exchange Process




The Cation Exchange Mechanism




Examples of Column Chemistries

Latex agglomerated Anion Exchanger




- ◆ Cross-Linked EVB-DVB core
- ◆ 5-10 μm diameter particles
- ◆ NanoBead™ anion-exchange layer

Hyper-Branched Surface Modified Anion Exchanger



Macro porous high capacity Anion Exchanger



Columns

Three factors that will affect resolution:

- ✓ Efficiency
- ✓ Capacity
- ✓ Selectivity

Altering Selectivity in Ion Chromatography

- Mobile phase parameters

- Eluent Ion
- Eluent concentration
- Eluent pH
- Organic additives
- Temperature

Mainly **retention** changes.


Little control on **selectivity**

- Stationary phase parameters

- Material
- Functional group
- Capacity


Better control of **selectivity**

Parameters Controlling Ion Exchange Retention



Effect of Stationary Phase Capacity

Number of ion exchange sites/weight equivalent of column packing.
The higher the capacity, the stronger the retention.



★ Most ion exchange columns have low capacity to improve detection and sensitivity

★ Complex samples and measurement of trace levels in the presence of high level matrix ions require higher capacities

Peaks:

1. Fluoride 2. Chloride 3. Nitrite
4. Phosphate 5. Bromide 6. Nitrate 7. Sulfate




Suppression Theory

- ✓ Suppression is defined as the neutralization by ion exchange of the salt of a weak acid or weak base eluent.
- ✓ Two results of suppression are:
 - ✓ Decreases background conductivity
 - ✓ Increases analyte response

Suppressed Conductivity Detection


- Batch Suppression
- Continuous Suppression

Batch Suppression with Conductivity Detection






Problem – Capacity and Dead Volume!

Batch-Mode Suppressor Configuration




Continuous Chemical Suppression




Continuous and low dead volume

Electronic Suppression




Ease of use with
lower backgrounds

Design of a Self Regenerating, Continuous Suppressor (SRS)



High Capacity with very small delay volumes


SRS® Ultra II AutoSuppression® Recycle Mode



Benefits of Continuous Suppression and the Important Link to Column Choice


- Larger Dynamic Working Range – High Sample load
- Allows the use of high capacity ion exchange columns
 - Use of high ionic strength mobile phases
- Allows the use of Gradient-Elution and higher resolution in combination with conductivity detection
- Coupling Ion Chromatography with other Detectors (e.g. MS, ICP-MS)

How is Conductivity Measured?



The conductivity of a solution is measured by applying an alternating voltage between two electrodes in a conductivity cell. At any instant in time, negatively charged anions migrate toward the positive electrode and positively charged cations migrate toward the negative electrode.

Isocratic Anion Separation



Column: IonPac® AG14A, AS14A, 5 μm , 3 x 150 mm

Eluent: 8.0 mM sodium carbonate
1.0 mM sodium bicarbonate

Flow Rate: 0.8 mL/min

Temperature: 30 ° C

Inj. Volume: 5 μL


Detection: Suppressed conductivity, ASRS® ULTRA II, 2 mm, AutoSuppression® recycle mode

Peaks:

1. Fluoride	5mg/L (ppm)
2. Acetate	20
3. Chloride	10
4. Nitrite	15
5. Bromide	25
6. Nitrate	25
7. Phosphate	40
8. Sulfate	30

Reproducibility of Common Anions on AS23 Column

Column: IonPac AS23, 4 × 250 mm
Eluent: 4.5 mM Na₂CO₃
0.8 mM NaHCO₃
Flow Rate: 1.0 mL/min
Inj. Volume: 25 µL
Detection: Suppressed conductivity




Peaks:

1. Fluoride	1.0 mg/L
2. Chloride	1.5
3. Nitrite	5.0
4. Bromide	5.0
5. Nitrate	5.0
6. Phosphate	7.5
7. Sulfate	7.5

Reproducibility of Common Cations in Drinking Water


Overlay of 300 Consecutive Runs




Column: IonPac CS12A, 4 × 250 mm
Eluent: 20 mM methanesulfonic acid
Flow Rate: 1.0 mL/min
Inj. Volume: 25 μL
Detection: Suppressed conductivity

Peak: 1. Sodium
2. Ammonium
3. Potassium
4. Magnesium
5. Calcium

Hydroxide Eluent Generation for Anion Analysis



A State of the art Ion Chromatography System





Improved Performance for Anion Gradient Separation

Overlay of 10 Consecutive Runs

Column: IonPac® AG11, AS11, 4 mm
Eluent: 0.5 to 25 mmol/L KOH
From: EG or 0.1 mol/L KOH


Flow Rate: 2 mL/min
Injection: 25 µL
Detection: Conductivity after ASRS®
Suppression, recycle mode



Peaks: 1. Fluoride 0.2 mg/L
2. Chloride 0.3
3. Nitrate 1.0
4. Sulfate 1.5
5. Phosphate 1.5



Purity and consistency of
Electronically Generated Eluents

High-Resolution Anion Analysis with Eluent Generation



Linearity of Sulfate with KOH Eluents

Waste sample
AS18 4 mm, 33 mM KOH



Fast, High Resolution Ion Chromatography

- Requirements
 - Fast, High Resolution Columns – Monoliths, small particle sizes
 - Low dead volume inert systems capable of higher pressures
 - Higher flow rates
 - Capillary systems

WHY?

- Faster more frequent analysis
- More information from the sample as resolution increases
- Always running and available for analysis
- Ease of use


Typical Conventional and Capillary IC Operating Parameters

	Conventional IC	Capillary IC
Column i.d.	4 mm	0.40 mm
Flow Rate	1.0 mL/min	10 μ L/min
Injection Loop	25 μ L	0.4 μ L
Suppressor Dead Volume	60 μ L	0.6 μ L
EG Current (50 mM KOH)	80.4 mA	0.804 mA
K ⁺ Consumption/Year	26.3 Moles (50 mM KOH)	0.263 Moles (50 mM KOH)
H ₂ O Consumption/Year	525.6 L	5.256 L

Separation of Common Anions using a Capillary AS19 Column – Can Capillary flow be reproducible


Peaks: 1. Fluoride 5. Chlorate
2. Bromate 6. Bromide
3. Chloride 7. Nitrate
4. Nitrite 8. Sulfate

Column: IonPac® AS19 column (0.4 mm × 25 cm)
Eluent Source: Capillary EGC-KOH cartridge
Eluent: 20 mM KOH
Flow Rate: 10 µL/min
Temperature: 30 ° C
Inj. Volume: 0.4 µL
Detection: Suppressed conductivity with ACES™ 300, AutoSuppression® recycle mode




Trace Level Anion Analysis on a Capillary Column with only 180µl injected

Column: IonPac® AS15 column (0.4 × 250 mm)
Eluent Source: EGC-KOH capillary cartridge
Eluent: 7 mM KOH (0 to 10 minutes), 7 to 32 mM KOH (10 to 16 minutes),
32 to 50 mM KOH (16 to 30 minutes), 50 to 65 mM
(30 to 33 minutes), 7 mM KOH (33 to 38 minutes)
Flow Rate: 12 µL/min
Temperature: 30 ° C
Detection: Suppressed conductivity, Anion Capillary Electrolytic Suppressor,
ACES™ 300, AutoSuppression recycle mode
Concentrator: IonSwift™ MAC-100 Concentrator (0.5 mm × 80 mm)
Injection Volume: 180 µL
Sample: A: Deionized water
B: Deionized water spiked with 0.5 µg/L Fluoride, 2.5 µg/L chloride,
nitrite, bromide, nitrate, and phosphate, 5.0 µg/L sulfate



Peak	A	B
1. Fluoride	0.018	0.48 µg/L
2. Chloride	0.12	2.49
3. Nitrite	0.042	2.53
4. Carbonate	—	—
5. Sulfate	0.075	4.72
6. Bromide	—	2.36
7. Nitrate	0.15	2.58
8. Phosphate	—	2.15


Fast IC Inorganic Capillary Anion Analysis at Different Flow Rates



Column: IonPac® AS18-Fast column (0.4 mm × 15 cm)
Eluent Source: Capillary EGC-KOH cartridge
Eluent: 33 mM KOH
Flow Rate: 10 to 20 µL/min
Temperature: 30 ° C
Detection: Suppressed conductivity with ACES™ 300, AutoSuppression® recycle mode
Inj. Volume: 0.4 µL

Peaks:	Retention Time at 20 µL/min (min)	Concentration (mg/L)
1. Fluoride	1.187	1.0
2. Chloride	1.623	1.5
3. Nitrite	1.937	5.0
4. Sulfate	2.197	7.5
5. Bromide	2.623	5.0
6. Nitrate	3.083	5.0
7. Phosphate	3.983	7.5

Fast IC: Inorganic Capillary Cation analysis at Different Flow Rates




Column: IonPac® CS16, 0.5 mm
Eluent Source: EGC-MSA Capillary Cartridge
Eluent: 30 mM MSA
Flow Rate: See Chromatogram
Inj. Volume: 0.4 μL
Column Temp.: 40 ° C
Detection: Suppressed Conductivity, Cation Capillary Electrolytic Suppressor, CCES™ 300 AutoSuppression recycle mode

Peaks:

1. Lithium	0.5	mg/L
2. Sodium	2.0	
3. Ammonium	2.5	
4. Potassium	5.0	
5. Magnesium	2.5	
6. Calcium	5.0	

Fast IC: Inorganic Capillary Cation analysis at Different Flow Rates




Column: IonSwift® MAX-200, 0.25 mm
Eluent Source: EGC-KOH capillary cartridge
Gradient: KOH gradient from 2 mM to 50 mM
Flow Rate: See Chromatogram
Inj. Volume: 0.4 µL
Column Temp.: 30 ° C
Detection: Suppressed Conductivity, Anion Capillary Electrolytic Suppressor, ACES™ 300, AutoSuppression recycle mode

- Peaks:
- | | |
|------------------|---------------------|
| 1. Fluoride | 11. Sulfate |
| 2. Acetate | 12. Oxalate |
| 3. Formate | 13. Tungstate |
| 4. Butyrate | 14. Phosphate |
| 5. Galacturonate | 15. Chromate |
| 6. Chloride | 16. Citrate |
| 7. Nitrite | 17. Isocitrate |
| 8. Bromide | 18. cis-Aconitate |
| 9. Nitrate | 19. trans-Aconitate |
| 10. Carbonate | |

Fermentation Broth Organic and Inorganic Anions


Using a monolithic 0.25 × 250 mm Column




Peaks:

1. Lactate	2 mg/L
2. Acetate	2
3. Propionate	2
4. Formate	2
5. Butyrate	2
6. Pyruvate	2
7. Isovalerate	2
8. Valerate	2
9. Chloride	2
10. Nitrate	2
11. Malate	2
12. Carbonate	—
13. Malonate	2
14. Maleate	2
15. Sulfate	2
16. Oxalate	2
17. Phosphate	4
18. Citrate	4
19. Isocitrate	4
20. Pyrophosphate	4

4µm Particle Size Fast Capillary Columns



Column: See chromatogram
Eluent: See chromatogram
Eluent Source: EGC-KOH capillary cartridge
Temperature: 30 ° C
Flow Rate: 30 µL/min
Inj. Volume: 0.4 µL
Detection: Suppressed conductivity,
Anion Capillary Electrolytic
Suppressor (ACES™ 300)
AutoSuppression® recycle mode



Peaks:	1.Fluoride	0.2	mg/L
	2.Chloride	0.5	
	3.Nitrite	1.0	
	4.Sulfate	1.0	
	5.Bromide	1.0	
	6.Nitrate	1.0	
	7.Phosphate	2.0	

4 μ m Particle Size High Resolution Columns

Column: IonPac AS11-HC-4 μ m HR,
capillary column

Eluent: 1-14 mM KOH in 16 min,
14-55 mM KOH in 24 min,

Eluent Source:

Flow Rate:

Inj. Volume:

Temperature:

Detection:

EGC KOH capillary cartridge


15 μ L/min

0.4 μ L

30 $^{\circ}$ C


Suppressed
conductivity, Anion
Capillary Electrolytic

ACESTM 300, 24 mA
AutoSuppression recycle
mode




- MS compatible Ion Exchange
 - Continuous Suppression IC system
- Electrolytic gradient generation
 - isocratic deionised water used from the pump
 - gradient generated pre-column electrolytically
- Electrolytic suppression
 - Exchange of post-column K^+ with H_3O^+ using electronic suppressor


IC-MS Schematics




IC separation of key metabolites that run poorly or not at all on HILIC






Ambiguity of identification due to fragmentation is common



ADP on 3 different stationary phases



Metabolomics Solution for charged intermediates – Capillary IC-MS



University Professor:
“IC separates what is in the
void volume of HPLC”

Summary

- **Wide range of IC Columns:**
 - Capacity, Selectivity, Specific applications
- **Suppression is the Key to Analytical Ion Chromatography**
 - Continuous Electrolytic Suppression allows the use of gradients and high capacity columns.
 - Robust, reproducible, quantitative and simple to use
- **Capillary IC a very powerful and easy to use analytical technique**
 - Always On/Always Ready, Ease of use, Fast with High Resolution
- **ICMS – Ultimate in High Performance, Extreme Flexibility and Information Rich for Method Development and Complex Samples**

Today's Webcast was brought to you by:



All questions from today's webcast will be answered in the CHROMacademy Forum

A screenshot of the LC|GC's CHROMacademy website. The header includes the logo and a search bar. A navigation menu contains links for Home, HPLC, GC, Sample Prep, Mass Spec, Ask the Expert, and Forum. The main content area features a 'FREE Membership' section with a list of benefits and a 'FREE Lite membership' button. To the right, there are circular icons for HPLC, GC, MS, and SPE. Below this is a section for 'CHROMacademy Essential Guide Webcast and Tutorials' featuring a webcast titled 'Translating GC Methods from Helium to Hydrogen Carrier Gas'. On the far right, a sidebar promotes the 'ESSENTIAL GUIDE WEBCAST' and includes a 'Got a problem? Ask our experts.' section with icons of three people and buttons for 'HPLC Troubleshooter', 'GC Troubleshooter', and 'The essential guides'.

www.chromacademy.com

Question and Answer

Type your question in the
“Submit Question”
box below your slide window.

Thank You

We Need Your Feedback!

Please complete the post-webcast evaluation form that will immediately appear in your window.